

The Development of a Leg Simulation Model to Assess the Fluid Handling Capabilities of Dressings



Michelle Delbono, Steve Foster, Risham Amjam, Systagenix, Gargrave, UK

ABSTRACT

Aim: To assess the fluid handling properties of dressings in a clinically relevant *in vitro* model.

Introduction: Wounds produce exudate which must be removed in order to prevent maceration and deterioration of the wound. Dressing performance is therefore often assessed by measuring total fluid handling capacity: a combination of absorbency and moisture vapour transmission rate. However, these tests may not correlate to how a dressing will perform clinically. We have therefore developed a simulated vertical leg model to assess how dressings manage wound fluid in a clinically-relevant situation.

Method: A vertical leg model was developed whereby simulated wound fluid containing a blue dye was continuously fed onto moist wound healing dressings (5.5–8 mL/24 hrs). Dressings were monitored daily for fluid absorption, retention and leakage.

Results: TIELLE® Max provided enhanced fluid handling capability compared with other moist wound healing dressings when tested in a clinically-relevant leg model. This was demonstrated by measuring the time taken for fluid to pool at the edge of the dressing: 4–6 days for TIELLE® Max compared with 24 hrs for the majority of other dressings tested. In addition, time until dressing leakage for TIELLE® Max was 5–7 days when a heavily exuding wound was simulated. Fluid was also distributed evenly throughout the TIELLE® Max dressing and retained, which will help to reduce the potential of wound maceration.

Conclusions: Results suggest that application of TIELLE® Max may result in fewer dressing changes than other moist wound healing dressings and will potentially reduce wound maceration.

RESULTS


Dressings at 24 hours

TIELLE® Max	ALLEVYN® Non-adhesive	MEPILEX®	BIATAIN® Soft Hold	BIATAIN® Non-adhesive	VERSIVA® XC Non Adhesive	ACTIVHEAL® Non-adhesive
TIELLE® Max and VERSIVA® did not show pooling of liquid at the edge of dressing until at least day 4 compared to all other dressings tested where pooling was observed after 24 hours.						


Dressings tested until failure

TIELLE® Max	ALLEVYN® Non-adhesive	MEPILEX®	BIATAIN® Soft Hold	BIATAIN® Non-adhesive	VERSIVA® XC Non Adhesive	ACTIVHEAL® Non-adhesive
Failure 5-7 days	Failure 2-6 days	Failure 3-4 days	Failure 2 days	Failure 1-4 days	Failure 5-6 days	Failure 1-7 days
<ul style="list-style-type: none"> •Best overall performance of all dressings tested •Longest time until breach of test fluid from the dressing •Minimal leakage onto test cylinder •Better distribution of fluid 	<ul style="list-style-type: none"> •Dressing performance inconsistent •Significant pooling of fluid at edge of dressing after 24 hours •Uneven distribution of fluid 	<ul style="list-style-type: none"> •Pooling of fluid at edge of dressing after 24 hours which then tracked upwards before failure 	<ul style="list-style-type: none"> •Short time before dressing failure •significant pooling of fluid at edge of dressing after 24 hours •Uneven distribution of fluid •Significant swelling at the back of the dressing 	<ul style="list-style-type: none"> •Dressing performance inconsistent •significant pooling of fluid at edge of dressing after 24 hours •Significant swelling at the back of the dressing 	<ul style="list-style-type: none"> •Long time until breach of test fluid from edge of dressing, however it became saturated and significant leakage of test fluid was observed on test cylinder. 	<ul style="list-style-type: none"> •Dressing performance inconsistent •Uneven distribution of fluid resulted in pooling of fluid at dressing edge after 24hours

Method



- The Leg Simulation model is designed to deliver a set quantity of test fluid into each dressing (5.5- 8ml/24hours)
- Dressings were secured to the cylinders with Sleek dressing tape
- The cylinders were kept upright for the length of the experiment
- Three replicates were tested for each dressing type
- Dressings were evaluated a minimum of every 24 hours until failure
- Dressing failure criteria was when the test fluid pooled outside the dressing, a large pooling of fluid within the dressing, or any other physical abnormalities.



TIELLE® Max provided enhanced fluid handling capability compared with other moist wound healing dressings tested on a clinically-relevant leg model

- TIELLE® Max consists of a hydro-polymer foam layer and a highly absorbent wicking layer
- The hydro-polymer layer absorbs wound exudate from the wound site and expands to fit the cavity of the wound
- The exudate is initially locked within the hydro-polymer, and is then transferred to a highly absorbent wicking layer enhancing the fluid handling capabilities

CONCLUSION

Results suggest that TIELLE® Max may result in fewer dressing changes when compared to the other moist wound healing dressings tested.

This is confirmed by the ability of TIELLE® Max to manage simulated wound fluid which was evenly distributed throughout the dressing during the experiment. The majority of other dressings tested in this study showed sub-optimal management of simulated wound fluid as demonstrated by pooling of the fluid at the dressing edge within 24 hours.

Wear time with TIELLE® Max was greatest, with time until dressing failure 5-7 days. This combined with minimum fluid remaining on the test cylinder, suggests optimal wound fluid management and reduced risk of wound maceration in a clinical setting.

* This product is a trademark of its owner

